



Carestream Health, Inc.
150 Verona Street
Rochester, NY 14608

RF Exposure Declaration

I, Dale B. Parks, Worldwide EHS manager of Carestream Health Inc., do hereby declare that the our product "KODAK 1500 Intraoral Camera" will operate with following manner, and will maintain RF exposure below 20mW (source-based time-averaged power) for all operation modes.


Signature: *Dale B. Parks*
Name: Dale B. Parks
Title: Manager World Wide EHS
Date: 08/25/2009



Boundary between inside/outside mouth for the camera



Power On/Off button



Position of antenna

Minimum distance of the antenna to person/mouth = 9 cm