

# Project file Module M21: 10.1 Instruction for use (IFU)

title	10.1 Instruction for use (IFU)			
responsible		creating system		
	Gerd Wotha	MS-Word		
author		document type		
	Gerd Wotha	doc		
filling		KM-tool		
	Design history file	tbd		
filling name and path				
I:\Techn Management\Projekte\Electronic_Accessory_Signature_System\DP-				
Project\05-Qualification\CE & FCC Akte\CE-Akte\10 Device Labelling\10.1 IfU incl accessory list\10.1 Instruction for use (IFU).doc				

### Released

name	date	signature
Gerd Wotha		

# 1. Instruction for use of Medical device

#### **1.1. Declaration of conformity**

A declaration of conformity according to R&TTE will be inserted to each instruction for use of the Medical device which uses the Module M21.

#### 1.2. General items

- M21 Reader Module is an integrated part within Draeger Medical devices. The antennas are permanently connected to the M21 via connectors.
- M21 Reader Module has an operational manual. The programming of the Draeger Medical device is based on this document. There is no user interface to the module.

#### 1.3. Mandatory text

- Changes or modifications not expressly approved by the party responsible for compliance could void the user's authority to operate the equipment
- The intended use of the M21 will be explained within each instruction for use of the Medical device which uses the Module M21.
- Following text has to be content of the device IFU: This device complies with Part 15 of the FCC Rules. Operation is subject to the following two conditions: (1) this device may not cause harmful interference, and (2) this device must accept any interference received, including interference that may cause undesired operation.