











Remote Control Artwork-Label-and -Location-VRT680135

Johnson and Johnson Surgical Version
Artwork file for
Part number VRT680135
FCC ID: VGEVRTREMOTE

GRAPHIC COLORS
NOT TEXTURES

P 83 - 16 C	
P 106 - 7 C	
P 118 - 15 C	
P 93 - 8 C	
P 173 - 10 C Lines and Solid Buttons	
Process Black C Main Background	
White C All Logos	
PMS 542C Bottle Blue	
COOL GREY 4C Bottle Grey	



KEYPAD EMOSS, 11PLS
.012H +/- .003
(EMOSS INCLUDE BUTTON BORDERS)

PART TEXTURE FINISHES
NOT COLORS







- BLACK Areas
Glossy Finish
- Grey Areas
Glossy Finish
- Blue Area
SPI TEXTURE
- Yellow Area
Clear with Heavy
Optical Diffuser

Johnson and Johnson Surgical Version
Label for VRT680135

Die line and dashed cut line does not print


Z353683

REF (Part No.): VRT680135			
SN (Serial Number):			
 (YYYY-MM-DD):			
Product of USA IPX4			XXXX
 R XXX-XXXXXX			Johnson & Johnson Surgical Vision, Inc. 1700 E. St. Andrew Place Santa Ana, CA 90705 USA


FCC ID: VGEVRTREMOTE IC ID: XXXXX-XXXXXXXX VERITAS™ Remote Control
 Remote Control Transmitter VRT680135 Patents: www.jnjvisionpro.com/patents

This device complies with part 15 of the FCC Rules. Operation is subject to the following two conditions:
 (1) This device may not cause harmful interference, and (2) this device must accept any interference received, including interference that may cause undesired operation.

Z353683 Rev. C 0320



Z353684



Johnson and Johnson Surgical Version
Label locations for VRT680135

