

4

3

2

1

NOTE:

1. USE ZEBRA PRINTER #TLP 2844 WITH 23-091-003 LABEL, DATA, ECOREX NS.BTW REV C DATA FILE

REVISIONS				
DCN	REV	DESCRIPTION	DATE	APPROVED
	C	Updated Description removed 2.0 Added RFID		

D

D

C

C

B

B

A


A

**EcoRex Collection Station**

Mfr:  
Vesta Medical, LLC  
Irvine, CA 92614

Cat. #  
Ser. #

~115 VAC, 1A  
60Hz  
FUSE: 250V 1A SLO-BLO 3AG



Read instructions  
for use before  
operating device

Contains FCC ID: XC7ECOREXRFID  
This device complies with part 15 of the FCC rules. Operation is subject to the following two conditions: (1) this device may not cause harmful interference and (2) this device must accept any interference that may cause undesirable operation.

ITEM	QTY.	PART NO.	DESCRIPTION
1	1	33-001-001	Label, Silver Polyester

DRAWN BY JORGE SERNA		DATE 9/17/08		<b>VESTA MEDICAL</b> IRVINE, CA (949) 660-8648 DESCRIPTION LABEL, DATA, RFID ECOREX			
CHECKED BY		DATE					
ENG APPR		DATE					
		SIZE B		DWG NO 23-091-003		REV C	
SCALE NONE			SHEET 1 OF 1				

4


3

2

1

NOTES UNLESS OTHERWISE SPECIFIED

1. MATERIAL: AVERY P/N 5267 OR EQUIVALENT
2. PRINT LABEL SHEET USING WORD TEMPLATE 33-100-001 REV B



UNLESS OTHERWISE SPECIFIED		NAME	DATE	<b>VESTA MEDICAL</b> IRVINE, CA (949) 660-8648
DIMENSIONS ARE IN INCHES TOLERANCES ANGULAR: ±1° TWO PLACE DECIMAL: ±.01 THREE PLACE DECIMAL: ±.005 BREAK ALL EDGES: .005 MAX SURFACE FINISH: 63		DRAWN	AH	
INTERPRET DRAWING PER ANSI Y14.5 1994		CHECKED		<b>LABEL, FCC, RFID IO PCB</b>
		ENG APPR		
MATERIAL	33-001-006	MFG APPR		
FINISH	NONE	QA		
DO NOT SCALE DRAWING		<b>CONFIDENTIAL INFORMATION</b> THE INFORMATION CONTAINED IN THIS DRAWING IS THE SOLE PROPERTY OF VESTA MEDICAL. ANY REPRODUCTION IN PART OR AS A WHOLE WITHOUT THE WRITTEN PERMISSION OF VESTA MEDICAL IS PROHIBITED.		SIZE <b>A</b>
		DWG NO <b>33-100-001</b>		REV <b>B</b>
		SCALE 2:1		SHEET 1 OF 1


B	ADDED FCC TWO PART STATEMENT			
REV	DESCRIPTION	DCN	DATE	APPR

FCC ID: XC7ECOREXRFID  
This device complies with Part 15 of the FCC Rules. Operation is subject to the following two conditions:  
1. This device may not cause harmful interference; and  
2. This device must accept any interference received, including interference that may cause undesired operation

I0brd\_V07  
24-065-001 x4

12-121-001

8







Data Label location