


FCC LABEL ID AND LOCATION

The Load Sensor is contained within a sterile package until immediately before use. The Load Sensor is activated while still in the sterile package, removed from the sterile package for actual use within a single surgery, and then both Load Sensor and packaging are discarded. Once activated, the Load Sensor has enough battery life for approximately 40 minutes of transmission, enough for one surgical use. The FCC ID is indicated on the label attached to the sterile packaging.





ORTHOSENSOR

1560 SAWGRASS PARKWAY, 4TH FLOOR
SUNRISE, FL 33323

