

DESCRIPTION

1. This document specifies the US FCC Part 15 Label for Bluetooth sensors to be affixed to the outside of the product box. Use this label in the following BOMs:
 - a. BOM-1076 VERASENSE for Stryker Triathlon GEN 4 (no int shims)
 - b. BOM-1071 VERASENSE for Stryker Triathlon GEN 4
 - c. BOM-1058 VERASENSE for Exactech Equinox Bill of Materials Naglreiter for US
 - d. BOM-1057 Stryker Mako Sensor Card BOM
 - e. BOM-1047 VERASENSE for Zimmer Biomet Persona (DreamTech)
 - f. BOM-1014 VERASENSE for Zimmer Biomet Persona (Creation Technologies)
 - g. BOM-1078 VERASENSE for Zimmer Biomet Persona with Alignment BOM

REFERENCE FIGURE

1. Label example shown below:

This device complies with part 15 of the FCC Rules. Operation is subject to the following two conditions: (1) This device may not cause harmful interference, and (2) this device must accept any interference received, including interference that may cause undesired operation.

LB-5248-0

NOTES

1. Print on OL2083 - 1.875" x .9375" Blank Label Template.
2. Print text as shown using black ink and Times New Roman font height of 6 point. Ensure that label document number and revision number are in the bottom right-hand corner.
3. Record all required information on DHR.

REVISION HISTORY

Rev	Description
0	Initial Release.