



# External Photos

**EUT Name:** VERASENSE for Zimmer Biomet Persona

**EUT Model:** CR C-D/3-9 Left, ZBH-PSNCRCD39-L, CR G-H/7-12 Left,  
ZBH-PSNCRGH712-L

**FCC ID:** XNL-ORTHOSNSR7

FCC Title 47, Part 15C, RSS-210 Issue 9, ANSI C63.10:2013

*Prepared for:*

Joshua Gardner  
OrthoSensor  
1855 Griffin Rd, St A-310  
Dania Beach, FL 33044 USA  
Tel: (954) 790-6744  
Fax: (954) 790-6744


*Prepared by:*

TUV Rheinland of North America  
762 Park Avenue  
Youngsville, NC 27596  
Tel: (919) 554-3668  
Fax: (919) 554-3542  
<http://www.tuv.com/>

*Report/Issue Date:* 8 June 2018

*Report Number:* Supplement to 31753637.001 - External Photos

## External Photos:

<h1>Memo</h1> <p><b>To:</b> FCC Validation <b>From:</b> Product Development <b>Date:</b> April 23, 2018 <b>Re:</b> Models for the VERASENSE for Zimmer Biomet Persona</p> <hr/> <p><b>Description:</b></p> <p>VERASENSE is indicated for any medical condition in which primary or revision Total Knee Arthroplasty (TKA) would be indicated. For use as a tool for adjustment of the femoral knee implant to reduce instability from flexion gap asymmetry. The VERASENSE is sterile, for single patient use.</p> <p><b>Sizes:</b></p> <p>Three sizes exist a Small, Medium, and Large, denoted as CD, EF, and GH respectively.</p> <p><b>Symmetry:</b></p> <p>Each size is offered for a left or right knee orientation denoted as "L" and "R" respectively.</p> <p>The internal electrical components of all models are identical. The size and symmetry are the only differences, a total of 6 models exist which belong to the VERASENSE for Zimmer Biomet Persona, see images below for examples:</p> <ul style="list-style-type: none"><li>• VERASENSE for Zimmer Biomet Persona CR C-D/3-9 Left-- ZBH-PSNCRCD39-L</li><li>• VERASENSE for Zimmer Biomet Persona CR C-D/3-9 Right-- ZBH-PSNCRCD39-R</li><li>• VERASENSE for Zimmer Biomet Persona CR E-F/3-11 Left-- ZBH-PSNCREF311-L</li><li>• VERASENSE for Zimmer Biomet Persona CR E-F/3-11 Right-- ZBH-PSNCREF311-R</li><li>• VERASENSE for Zimmer Biomet Persona CR G-H/7-12 Left-- ZBH-PSNCRGH712-L</li><li>• VERASENSE for Zimmer Biomet Persona CR G-H/7-12 Right --ZBH-PSNCRGH712-R</li></ul>	
--	--


**Image 1:** CR C-D/3-9 Left-- ZBH-PSNCRCD39-L, CR E-F/3-11 Left-- ZBH-PSNCREF311-L, and CR G-H/7-12 Left --ZBH-PSNCRGH712-L



**Image 2:** CR C-D/3-9 Left-- ZBH-PSNCRCD39-L



**Image 3:** CR E-F/3-11 Left-- ZBH-PSNCREF311-L



**Image 4:** CR G-H/7-12 Left --ZBH-PSNCRGH712-L





Figure 1: Front view of model CR G-H/7-12 LEFT, ZBH-PSNCRGH712-L and model CR C-D/3-9 LEFT, ZBH-PSNCRCD39-L.

Both models were investigated.

The Model CR G-H/7-12 LEFT, ZBH-PSNCRGH712-L was selected for testing

See Clause 4 of the test report




Figure 2: Rear view of model CR G-H/7-12 LEFT, ZBH-PSNCRGH712-L and model CR C-D/3-9 LEFT, ZBH-PSNCRCD39-L.

Both models were investigated.

The Model CR G-H/7-12 LEFT, ZBH-PSNCRGH712-L was selected for testing

See Clause 4 of the test report