



Intel® Clane User Manual

DRAFT

Home Deployment
August 31, 2009

Revision Number: 0.3

1. Important Safety Information and Warnings

| | |
|--|---|
| FCC Part 15 Compliance Statement..... | 3 |
| Information for the user..... | 3 |
| Explosive Device Proximity Warning | 4 |
| Safety Precautions..... | 4 |

2. Getting Started

| | |
|--------------------------|---|
| Introduction | 5 |
| Starting a Session | 5 |
| Ending a Session | 5 |

1. Important Safety Information and Warnings

FCC Part 15 Compliance Statement

This equipment has been tested and found to comply with the limits for a Class B digital device, pursuant to part 15 of the FCC Rules. These limits are designed to provide reasonable protection against harmful interference in a residential installation. This equipment generates, uses and can radiate radio frequency energy and, if not installed and used in accordance with the instructions, may cause harmful interference to radio communications. However, there is no guarantee that interference will not occur in a particular installation. If this equipment does cause harmful interference to radio or television reception, which can be determined by turning the equipment off and on, the user is encouraged to try to correct the interference by one or more of the following measures:

- Reorient or relocate the receiving antenna.
- Increase the separation between the equipment and receiver.
- Connect the equipment into an outlet on a circuit different from that to which the receiver is connected.
- Consult the dealer or an experienced radio/TV technician for help.
- Any changes or modifications not expressly approved by the grantee of this device could void the user's authority to operate the equipment. The customer is responsible for ensuring compliance of the modified product.
- This device complies with part 15 of the FCC Rules. Operation is subject to the following two conditions: (1) This device may not cause harmful interference, and (2) this device must accept any interference received, including interference that may cause undesired operation.

Information for the user

"Intel Corporation and Clane"

Safety Notices

The FCC with its action in ET Docket 96-8 has adopted a safety standard for human exposure to radio frequency (RF) electromagnetic energy emitted by FCC certified equipment. "Intel Corporation and Clane" meets the Human Exposure limits found in OET Bulletin 65, 2001, and ANSI/IEEE C95.1, 1992. Proper operation of this radio according to the instructions found in this manual will result in exposure substantially below the FCC's recommended limits. The following safety precautions should be observed:

- Do not touch or move antenna while the unit is transmitting or receiving.
- Do not hold any component containing the radio such that the antenna is very close or touching any exposed parts of the body, especially the face or eyes, while transmitting.

Use in specific environments:

The use of wireless devices in hazardous locations is limited by the constraints posed by the safety directors of such environments.

- The use of wireless devices in hospitals is restricted to the limits set forth by each hospital.

Antenna use:

- In order to comply with FCC RF exposure limits, low gain integrated antennas should be located at a minimum distance of 20cm (7.87 inches) or more from the body of all persons.

Explosive Device Proximity Warning

- Warning: Do not operate a portable transmitter (such as a wireless network device) in an explosive environment unless the device has been modified to be qualified for such use.

Safety Precautions

This product is professionally installed and not to be moved by users. Please read, follow, and keep these instructions and observe all warnings. **WARNING :**

- Do not open the Clane sensors. There are no user serviceable parts inside.
- Do not spill liquid on the Clane sensors.
- Do not operate the Clane sensors near water.
- Do not wipe the Clane sensors with a wet cloth.

2. Getting Started

Introduction

The Clane system transforms eldercare with advanced sensor technology that learns the daily living patterns of senior residents. The care staff sends you a series of questions that facilitates the delivery of personalized care to help improve response times and identify potential problems before they become emergencies.

Starting a Session

When you first get the prompt on your touchscreen, simply touch the screen to start the session and follow the instructions on the screen.

- For any question if, you are not ready to respond, you may defer the session by pressing the "Ask Later" button on the screen and you will be prompted later.
- When ready, answer each question by pressing the appropriate button (Yes or No) on the screen.

Ending a Session

- After the session, leave the touchscreen on.

Note: Please inform your care staff of your travel schedule when you leave the house for more than one day.



©2009, Intel Corporation. All rights reserved. Intel and the Intel logo are trademarks or registered trademarks of Intel Corporation or its subsidiaries in the United States and other countries.

*Other names and brands may be claimed as the property of others.