

Spinal Modulation  
1135 O'Brien Dr.  
Menlo Park, CA 94025

<http://www.spinalmodulation.com/>

Direct: 650 543-6829

Fax: 650 327-2336

15March2014      **MN20600-02 Tune-Up Procedure**

Spinal Modulation Patient Programmer Product MN20600-02.

The Basestation Board is trimmed at time of manufacturing which remains in effect the entire lifetime of the device. No tune-up procedure exists or is required.

April Pixley  
Senior Director, Hardware/Firmware Engineering  
Chief Systems Engineer

