kGOAL

Therapy Holdings, Inc. Boost L4-010 ID Label Information and Location (FCC)

The equipment "can only operate in conjunction with another associated product that has an electronic display" (it is controlled via an app on a smartphone and has no user controls on the device itself). Per FCC and IC guidance, e-labeling is valid for the product.

Please see below screen shot which shows where the app lists FCC ID and model number. The settings page, where the information is listed, is reachable within 1 click of the app's Home / main page.

 ✓ Search III 	? 4:4	2 PM	94% 🛃	2
	SETT	INGS		
Welco	me Flow			
язык			English	
App V	ersion		1.0	
Terms	of Service			
Privac	y Policy			
	D: 2A3MP-L4 I: L4-001	001		
Dashboard	Workout	Settings	P	

For software security, the FCC ID text is stored in the app (hard coded and cannot be modified or overwritten without publishing a new version of the app through the Apple App Store or Google Play Market). The app is exclusively for use with this Boost L4-010 product.

kGOAL

The below figure shows the physical temporary ID label on the product packaging per KDB 784748 D02 clause 4.5.

Designed in California. Made in China. US and international patents pending. The Boost and kGoal names, logos, and taglines are trademarked.



