




Details of: Digital Automatic Blood Wrist Pressure Monitor Internal view 1




Details of: Digital Automatic Blood Wrist Pressure Monitor Internal view 2



Details of: Digital Automatic Blood Wrist Pressure Monitor PCB view 1




Details of: Digital Automatic Blood Wrist Pressure Monitor PCB view 2







Details of: Digital Automatic Blood Wrist Pressure Monitor PCB view 3




Details of: Digital Automatic Blood Wrist Pressure Monitor PCB view 4



Details of: Digital Automatic Blood Wrist Pressure Monitor PCB view 5




Details of: Digital Automatic Blood Wrist Pressure Monitor PCB view 6







Details of: Digital Automatic Blood Wrist Pressure Monitor PCB view 7




Details of: Digital Automatic Blood Wrist Pressure Monitor PCB view 8



Details of: Digital Automatic Blood Wrist Pressure Monitor pump view 1



Details of: Digital Automatic Blood Wrist Pressure Monitor motor view 1





Details of: Digital Automatic Blood Wrist Pressure Monitor motor view 2

