




# Health & Life Co. Ltd.

FCC ID: 2ABTAHNL85A1Z

The FCC ID will be placed on the equipment as shown in the photograph below.

<b>Automatic Upper Arm Blood Pressure Monitor</b>	
Model No. : HL858A1-Z	
Power :	
1. AAA "LR03" (1.5V) Alkaline Battery x 4	
2. 5V 1A adapter(for US): <u>FranMar International, FRM06-S05-UU</u>	
(optional) 100-240V, AC, 50-60 Hz, 0.2 A	
5V 1.2A adapter(for EU) : <u>Zhongshan Baolijin Electronic Co., Ltd,</u>	
BLJ06L050120U-R (optional) 100-240V, AC, 50-60 Hz, 0.2 A	
-----	
Rated Range of Cuff Pressure : 0 ~ 300 mmHg	
Rated Range of Determination : 40 ~ 280 mmHg	
FCC ID: 2ABTAHNL85A1Z	
<b>SN</b> (XX)XXXXXXXXXXXXXXXX (XX)YYMMXXXXXX	Manufacturer : <b>HEALTH &amp; LIFE CO., LTD.</b> 9F, No. 186, Jian Yi Road, <u>Zhonghe District, New Taipei City, Taiwan</u>
   <b>IP22</b>	Made in China
This device complies with Part 15 of the FCC Rules. Operation is subject to the following two conditions: (1) this device may not cause harmful interference, and (2) this device must accept any interference received, including interference that may cause undesired operation.	
Date of Manufacture : YYYY-MM	