

FCC & ISED Labelling Submission Air 11

ResMed Confidential

Proprietary Information, not to be reproduced or made available to third parties without prior consent from ResMed and not to be used in any unauthorised way. © 2021 ResMed.

SPEC NO.	eLB1392-074
REV	3.0
PAGE	1 of 3

Melanie Webster 18 February 2021

Common labels/moulded in artwork



BRAND


ResMed branding chrome or colour foiled on <u>front</u> of product



"BF APPLIED PART"
Body Fibrillation Icon
Moulded into plastic case
using an insert
on rear of device,
adjacent to the device outlet




"R390-411 Label - MR Unsafe"
White pre-printed label stuck next to power connector on rear of assembly



"R390-413 Label - FCC & IC NIR"

White pre-printed label stuck on bottom of assembly





A390-7374 Bottom Case Assembly Laser-Etched on bottom of chassis



FCC & ISED Labelling Submission Air 11

ResMed Confidential

Proprietary Information, not to be reproduced or made available to third parties without prior consent from ResMed and not to be used in any unauthorised way. © 2021 ResMed.

SPEC NO.	eLB1392-074
REV	3.0
PAGE	2 OF 3

Melanie Webster 18 February 2021

Product Specific Labels for FCC













NOTE:

- The device ID is the base code, which can be packaged in different boxes of accessories
- CO means it is packed with a humidifier.
- TRI means it is packed with a humidifier and heated tube.



FCC & ISED Labelling Submission Air 11

ResMed Confidential

Proprietary Information, not to be reproduced or made available to third parties without prior consent from ResMed and not to be used in any unauthorised way. © 2021 ResMed.


SPEC NO.	eLB1392-074
REV	3.0
PAGE	3 OF 3

Melanie Webster 18 February 2021

Product Specific Labels for ISED













NOTE:

- The device ID is the base code, which can be packaged in different boxes of accessories
- CO means it is packed with a humidifier.
- TRI means it is packed with a humidifier and heated tube.