

AIRSENSE 11 - US VARIANTS - SIMILARITIES AND DIFFERENCES FOR PACIFIC PROJECT

## *Title:* AIRSENSE 11- US VARIANTS – SIMILARITIES AND DIFFERENCES

Project/Product PACIFIC

ResMed Confidential Proprietary Information, not to be reproduced or made available to third parties without prior consent from ResMed and not to be used in any unauthorised way. © ResMed Ltd 2022

## 1. Purpose

This document discusses the similarities and differences between US variants of the AirSense 11 project applicable to the following device models:

Model Name / Product Code	Marketing Name
39517	AirSense 11 AutoSet USA
39518	Airsense 11 CPAP USA
39519	AirSense 11 Elite USA
39520	AirSense 11 AutoSet CAN
39521	AirSense 11 CPAP CAN
39522	AirSense 11 Elite CAN

## 2. Equivalence

All the device models mentioned in <u>Section 1</u> have the following <u>same</u> components/parts:

- Hardware
  - o Cellular
  - o Bluetooth
  - Flow Generator
- Software
  - o Cellular
  - o Bluetooth
  - Flow Generator
  - Device Colour (paint)

## 3. Differences

Different models are derived for <u>marketing purposes</u>. The difference between Elite, AutoSet, and CPAP models are the therapy modes they support, whereas different variants of AutoSet, Elite, and CPAP are based on the customer centric features supported by the device. These settings are configured at manufacturing via a configuration file and is independent of the cellular, Bluetooth and FG software.

In summary, the differences between the device models mentioned in <u>Section 1</u> are:

- Therapy modes
- Data settings
- UI language
- Name plaque

**Company Confidential**