

The statement of model difference

We, Shenzhen Viatom Technology Co., Ltd. do hereby certify that, there are 4 models of Pulse Oximeter (PO5,PO5-1,PO5-2,PO5-3) in this registration. The only difference between the models is the color of main unit shell, the appearance structure, circuit structure, performance index and other functions are the same. We choose PO5 to be the typical test sample, which tested is the most complex and advanced model and represents the safety and efficiency of the models listed below.

The differences between models are as follows:

No.	Model	Appearance color	Intended users	Measuring part	Wavelength	Prompt tone
1	PO5	Blue	Adults, children and newborns	Finger/Foot	Visible red Light: 660nm Infrared light : 940nm	Buzzer
2	PO5-1	Pink				
3	PO5-2	White				
4	PO5-3	Green				

Regards,

Signature: 

printed name: Bing Wang