

PHOTOGRAPHS OF THE TEST SET-UP

Photo 1: Measurement System DASY5



Photo 2: Front View



Photo 3: Rear View



Photo 4: Left hand touched




Photo 5: Left hand tilted 15°




Photo 6: Right hand touched




Photo 7: Right hand tilted 15°




Photo 8: Front Side 15mm



Photo 9: Back Side 15mm



Photo 10: Front Side 10mm



Photo 11: Back Side 10mm



Photo 12: Left Side 10mm



Photo 13: Right Side 10mm	Photo 14: Top Side 10mm
	
Photo 15: Bottom Side 10mm	N/A
	N/A

Photograph: Liquid depth









Photo 16: HSL850 Depth (15.1cm)	Photo 17: HSL1900 Depth (15cm)
	
Photo 18: HSL2450 Depth (15.0cm)	Photo 19: Body 850 Depth (15.0cm)
	
Photo 20: Body1900 Depth (15.0cm)	Photo 21: Body2450 Depth (15.0cm)
	

Photo 22: Head2600 Depth(15.0cm)	Photo 23: Body2600 Depth(15.0cm)
	
Photo 22: Head1750 Depth(15.0cm)	Photo 23: Body1750 Depth(15.0cm)
