

LT5011 Box Sleeve

145mm

75mm

145mm

75mm

160mm

ThermaCare®



THERARELIEF DEVICE

Pain Therapy

Wearable Pain Relief TENS Technology

The Wireless Pain Relieve Device LT5011C is designed to be used for temporary relief of pain associated with sore and aching muscles in the shoulder, waist, back, upper extremities (arm) and lower extremities (leg) due to strain from exercise or normal household work activities. And to be used for the symptomatic relief and management of chronic, intractable pain and relief of pain associated with arthritis.


For all contraindications, warnings, precautions, and adverse side effects and proper use, please refer to the TheraRelief device user manual.



LOW BACK PAIN



SHOULDER PAIN



ELBOW PAIN

Includes:
TheraRelief Unit with rechargeable lithium ion battery
Micro-USB cable
Remote control with 2 AAA Batteries included
3 Adhesive electrode pads

This device complies with Part 15 of the FCC Rules. Operation is subject to the following two conditions: (1) this device may not cause harmful interference, and (2) this device must accept any interference received, including interference that may cause undesired operation.

Made in China

ThermaCare®

THERARELIEF DEVICE

Pain Therapy



Wearable Pain Relief TENS Technology

Relief of muscle, joint and back pain

Chronic pain relief

This device complies with Part 15 of the FCC Rules. Operation is subject to the following two conditions:


(1) this device may not cause harmful interference, and (2) this device must accept any interference received, including interference that may cause undesired operation.



KNEE PAIN



SHOULDER PAIN



LEG PAIN

Manufactured for PFIZER INC.
Add: PO BOX 34600 Bartlett
TN 38184-0600 United States
www.pfizer.com



UDI

LOT

