



FCC ID: 2AS8ABPMF001

This device complies with Part 15 of the FCC Rules. Operation is subject to the following two conditions: (1) this device may not cause harmful interference, and (2) this device must accept any interference received, including interference that may cause undesired operation.

size: 20*15mm

Product Name: Blood Pressure Monitor Model No.: F1701TRated input: 5.0V == 1.0A



Battery: DC 4.5V(3×1.5V AA batteries)

Shenzhen Jamr Technology CO., Ltd A101-301,D101-201, Jamr Science & Technology Park, No. 2 Guiyuan Road, Guixiang Community,Guanlan Street, Longhua District, 518100 Shenzhen, PEOPLE'S REPUBLIC OF CHINA

Shanghai International Holding Corp. GmbH (Europe) |P21 | Eiffestrasse 80, 20537 Hamburg, Germany

202104210F1701T0001



size: 39.8*15mm