


size:33*19mm



For in vitro diagnostic use and self-testing only. The meter is intended to be used for the quantitative measurement of glucose in fresh capillary whole blood with the Sejoy BS-362 Blood Glucose Test Strips.

SN XXXXXXXXXXX Manufacture date:xxxxxxx

FCC ID: 2AWPL0001
Manufactured by Hangzhou Sejoy Electronics & Instruments CO.,LTD.
Area C., Building 2. No. 365, Wuzhou Road, Yuhang Economic
Development Zone Hungzhou City 311100 Zhejiang China

Battery: ___ DC3.0V 1×CR2032











