

dēxter®

Dexter L6 - Patient Cart 2



UDI



REF 104

SN 12345

MD



230kg

(01)07640184840011  
(11)230801  
(21)12345

CE 1639



2023-08-01



700 VA  
100 - 240 V~  
50/60 Hz



T 10A H  
250 V~



802263



Distalmotion SA  
Route de la Corniche 3B  
1066 Epalinges  
Switzerland

EC REP

Qarad EC-REP BV  
Pas 257  
2440 Geel  
Belgium

FCCID: 2BCEADM-L6001

This device complies with Part 15 of the FCC Rules.  
Operation is subject to the following two conditions:

(1) this device may not cause harmful interference, and  
(2) this device must accept any interference received,  
including interference that may cause undesired  
operation.

L-112-01