

APPENDIX D: SAR TISSUE SPECIFICATIONS

Measurement Procedure for Tissue verification:

- 1) The network analyzer and probe system was configured and calibrated.
- 2) The probe was immersed in the tissue. The tissue was placed in a nonmetallic container. Trapped air bubbles beneath the flange were minimized by placing the probe at a slight angle.
- 3) The complex admittance with respect to the probe aperture was measured
- 4) The complex relative permittivity ε can be calculated from the below equation (Pournaropoulos and Misra):

$$Y = \frac{j2\omega\varepsilon_{r}\varepsilon_{0}}{[\ln(b/a)]^{2}} \int_{a}^{b} \int_{a}^{b} \int_{0}^{\pi} \cos\phi' \frac{\exp[-j\omega r(\mu_{0}\varepsilon_{r}'\varepsilon_{0})^{1/2}]}{r} d\phi' d\rho' d\rho$$

where Y is the admittance of the probe in contact with the sample, the primed and unprimed coordinates refer to source and observation points, respectively, $r^2 = \rho^2 + \rho'^2 - 2\rho\rho'\cos\phi'$, ω is the angular frequency, and $j = \sqrt{-1}$.

nts:	
Ethanediol	>1.0-4.9%
STOT RE 2, H373;	
Acute Tox. 4, H302	
Sodium petroleum sulfonate	< 2.9%
Eye Irrit. 2, H319	
•	
Hexylene Glycol / 2-Methyl-pentane-2,4-diol	< 2.9%
Skin Irrit. 2, H315; Eye Irrit. 2, H319	
Alkoxylated alcohol, > C ₁₆	< 2.0%
Aquatic Chronic 2, H411:	
ses refer to section 16.	
	Acute Tox. 4, H302 Sodium petroleum sulfonate Eye Irrit. 2, H319 Hexylene Glycol / 2-Methyl-pentane-2,4-diol Skin Irrit. 2, H315; Eye Irrit. 2, H319 Alkoxylated alcohol, > C ₁₆ Aquatic Chronic 2, H411; Skin Irrit. 2, H315; Eye Irrit. 2, H319

Figure D-1

Note: Liquid recipes are proprietary SPEAG. Since the composition is approximate to the actual liquids utilized, the manufacturer tissue-equivalent liquid data sheets are provided below.

FCC ID: BCGA2759	WIFI 6 GHZ RF EXPOSURE EVALUATION	Approved by: Technical Manager	
DUT Type: Tablet Device		APPENDIX D: Page 2 of 3	

-5.2

-0.8

-0.7 -0.5

0.6

0.8

1.4

1.6 2.0

2.1

2.0

1.6

1.2

Zeughausstrasse 43, 8004 Zurich, Switzerland Phone +41 44 245 9700, Fax +41 44 245 9779 www.speag.swiss, info@speag.swiss

Measurement Certificate / Material Test

Item Name Head Tissue Simulating Liquid (HBBL600-10000V6)

Product No. SL AAH U16 BC (Batch: 210629-3)

Manufacturer SPEAG

Measurement Method

TSL dielectric parameters measured using calibrated DAK probe.

Target Parameters

Target parameters as defined in the IEEE 1528 and IEC 62209 compliance standards.

Test Condition

Ambient Condition 22°C; 30% humidity

TSL Temperature 22°C Test Date 1-Jul-21 Operator WM

Additional Information
TSL Density TSL Heat-capacity

Result	s							Sec. 1							
	Measu	ured		Targe	t	Diff.to Tar	get [%]	15,0							
[MHz]	e'	e"	sigma	eps	sigma	∆-eps	∆-sigma	10.0					186		13
300	44.7	25.5	0.85	42.7	0.88	4.6	-3.6	% 5.0							
750	44.1	21.6	0.90	41.9	0.89	5.1	0.7								
00	44.0	20.6	0.92	41.7	0.90	5.6	2.5	Permittivity 0.0							
25	44.0	20.2	0.93	41.6	0.91	5.8	2.6	-5.0							
35	44.0	20.0	0.93	41.5	0.91	5.9	2.0	3-10.0 -15.0	- 1				4 10		_
50	43.9	19.8	0.93	41.5	0.92	5.8	1.5		00 150	00 0500	2500 45	00 5500 0	500 7500	0500.0	
0	43.8	19.0	0.95	41.5	0.97	5.5	-2.1	,	000 150	JU 2500		icy MHz	500 7500	8500 9	50
00	42.8	15.1	1.18	40.6	1.18	5.4	0.0	15.0							Ξ
50	42.7	14.9	1.20	40.5	1.20	5.4	0.0	10.0	90			MURRI			H
00	42.4	14.4	1.28	40.3	1.28	5.2	-0.3	8		A				348	Ī
25	42.4	14.3	1.30	40.3	1.30	5.3	0.1	Conductivity 200	A	11					
10	42.4	14.3	1.31	40.3	1.31	5.3	0.3	ducti	10	1					
0	42.3	14.3	1.31	40.2	1.31	5.1	-0.2	00-5.0							
0	42.3	14.2	1.34	40.2	1.34	5.3	-0.2	910.0 2015.0			tion i	4	1		Ī
0	42.2	14.1	1.37	40.1	1.37	5.3	-0.1		00 150	0 2500 :	3500 450	00 5500 6	500 7500	8500 9	50
)	42.1	14.0	1.40	40.0	1.40	5.3	0.0					ncy MHz			_
)	42.1	13.9	1.41	40.0	1.40	5.3	0.7	3500	39.4	14.2	2.77	37.9	2.91	3.8	_
3	42.1	13.9	1.42	40.0	1.40	5.3	1.4	3700	39.0	14.4	2.96	37.7	3.12	3.6	
	42.0	13.9	1.43	40.0	1.40	5.0	2.1	5200	36.4	16.0	4.62	36.0	4.66	1.2	
	42.0	13.8	1.46	40.0	1.40	5.0	4.3	5250	36.3	16.0	4.68	35.9	4.71	1.1	
	41.9	13.8	1.49	40.0	1.40	4.7	6.4	5300	36.2	16.1	4.73	35.9	4.76	1.0	
)	41.8	13.7	1.53	40.0	1.40	4.5	9.3	5500	35.9	16.2	4.96	35.6	4.96	0.6	
	41.8	13.7	1.56	39.9	1.44	4.7	8.0	5600	35.7	16.3	5.08	35.5	5.07	0.4	
	41.7	13.7	1.59	39.8	1.49	4.7	6.8	5700	35.5	16.4	5.20	35.4	5.17	0.2	
	41.6	13.6	1.63	39.7	1.53	4.7	6.3	5800	35.3	16.5	5.31	35.3	5.27	0.0	
)	41.6	13.6	1.67	39.6	1.58	4.9	5.8	6000	34.9	16.6	5.55	35.1	5.48	-0.4	
0	41.5	13.6	1.70	39.6	1.62	4.9	4.8	6500	34.0	17.1	6.17	34.5	6.07	-1.3	
)	41.4	13.6	1.74	39.5	1.67	4.9	4.4	7000	33.1	17.4	6.78	33.9	6.65	-2.2	
0	41.3	13.6	1.78	39.4	1.71	4.9	4.0	7500	32.3	17.7	7.40	33.3	7.24	-3.1	
0	41.3	13.6	1.82	39.3	1.76	5.1	3.7	8000	31.4	18.0	8.01	32.7	7.84	-4.1	
0	41.2	13.6	1.86	39.2	1.80	5.1	3.3	8500	30.5	18.2	8.62	32.1	8.45	-5.0	
)	41.1	13.6	1.90	39.1	1.85	5.0	2.5	9000	29.7	18.4	9.22	31.5	9.08	-5.9	
)	41.0	13.7	1.94	39.1	1.91	4.9	1.6	9500	28.9	18.6	9.82	31.0	9.71	-6.7	
0	41.0	13.7	1.98	39.0	1.96	5.1	0.8	10000	28.1	18.7	10.42	30.4	10.36	-7.6	

Figure D-2 600 - 10000 MHz Head Tissue Equivalent Matter

FCC ID: BCGA2759	WIFI 6 GHZ RF EXPOSURE EVALUATION	Approved by: Technical Manager	
DUT Type: Tablet Device		APPENDIX D: Page 3 of 3	