Nedap SmartTag-B Neck series

Label Information

Version	0.1
Date	02-12-2022, rev. 22-06-2023
Author	K. Damhuis
Checked by	A.P. Haytema

Model differences

The models are identical w.r.t. pcb lay-out, component placement and transmitting circuitry. Difference is made by software features and optional identification LF ferrite coil.

SmartTag-B Neck series, model breakdown

		Features				
Model	Туре	Identification	Fertility	Eating	Rumination	
SmartTag-B Neck	FE4		Х	Х		
	IFE4	Х	Х	Х		
	FER4		Х	Х	Х	
	IFER4	Х	Х	Х	Х	

Model **SmartTag-B Neck IFER4** is considered as representative for this document since it contains the maximum features.

Front of the Tag





page 1

Location of the label - On backside of the assembled tag.



Sample of label (all models)



Typedesignation (1 example) - laser imprinted on front



