

Return address: PO Box 6, 7140 AA Groenlo, The Netherlands

ATCB
Attn.: Mrs. M. Bosley
Certification Department
6731 Whittier Avenue, Suite C110
McLean, Virginia 22101
USA

Groenlo, 26 October 2011

Subject: Attestation of Similarity FCC ID: CGDRTLACT

Dear Mrs. Bosley,

The device with Model: LACTIVATOR REALTIME as mentioned in report "12_11060901.fcc01_Lact RT_Testreport_FCC-IC CGDRTLACT.pdf" is a transmitter operating on 433 MHz to detect activity of cows.

Model: REALTIME LACTIVATOR NEDAP is delivered by Nedap directly to its dealers and/or end users. Nedap also delivers 5 more models which are electronically as well as mechanically the same as the Model: LACTIVATOR REALTIME NEDAP only the color differs and it is delivered to other companies as mentioned in the list below who in turn deliver it to their dealers and/or end users, so another name appears on the product. Therefore more following Model numbers are introduced.

- | | |
|---|-------------------------------|
| 1. Model: LACTIVATOR REALTIME NEDAP | is yellow with Nedap |
| 2. Model: REALTIME LACTIVATOR LELY | is red with LELY |
| 3. Model: REALTIME LACTIVATOR GEA | is green with GEA |
| 4. Model: REALTIME LACTIVATOR BOUMATIC | is dark blue with BouMatic |
| 5. Model: REALTIME LACTIVATOR INTERPULS | is medium blue with InterPuls |
| 6. Model: REALTIME LACTIVATOR CRV | is light blue CRV |

All the products will of course be labeled with FCC ID: CGDRTLACT. The Model number will appear on the box containing 10 products each.

So the labeling of the product will be:

Model: LACTIVATOR REALTIME

FCC ID: CGDRTLACT


The product can be recognized as shown in the picture below by color and by Brand name and combined as mentioned in the table above.



I thereby wish to request whether you can add these models to the list of approved systems.

I hope you can honor our request

Best regards,
N.V. Nederlandsche Apparatenfabriek "Nedap"


Jacques Hulshof
Approbation Management