Each units TX and RX functionality is tested and judged with a pass/fail criteria during production. Units meeting the criteria continue through production and others are returned for diagnostics.

Rogers Labs, Inc. 4405 West 259<sup>th</sup> Terrace Louisburg, KS 66053

Revision 1

Dynastream Innovations, Inc.

Model: HRM2A Test #: 090614 SN: 33393 Phone/Fax: (913) 837-3214 Test to: CFR47 Parts 2 and 15.249, RSS-210 TunPro HRM2A

FCC ID#: O6RHMR2 IC: 3797A-HRM2 Page 1 of 1

Date: June 23, 2009