

DF4/IS4

Acticor 7 HF-T QP

REF 429522

SN 12345678

PID 89

IS-1
IS4
DF4-LLHH

2020-04-30

113 °F




1010403547915667117120043012112345678





Acticor 7 HF-T QP

ProMRI

PID 89



402-405 MHz
TP2



441732

FCC-ID: QRI-ICD4200

Distributor:
BIOTRONIK, Inc.
6024 Jean Road
Lake Oswego, OR 97036-5369
(800) 547-0394 (24-hour)
(503) 635-9936 (Fax)
sales@biotronik.com

Manufacturer:
BIOTRONIK SE & Co. KG
Woermannkehre 1
12359 Berlin - Germany
Tel. +49 (0) 30 68905-0
Fax +49 (0) 30 6852804
sales@biotronik.com

BIOTRONIK

Acticor 7 HF-T QP DF4/IS4

REF 429522 SN 12345678 PID 89




BIOTRONIK

Acticor 7 HF-T QP DF4/IS4
SN 12345678
PID 89
BIOTRONIK

Acticor 7 HF-T QP DF4/IS4
SN 12345678
PID 89
BIOTRONIK

Acticor 7 HF-T QP

IS-1
IS4
DF4-LLHH



DF4/IS4
ProMRI

REF 429522

2020-04-30


SN 12345678

PID 89

STERILEEO



Rx only



BIOTRONIK SE & Co. KG
Woermannkehre 1
12359 Berlin - Germany

BIOTRONIK



1010403547915667117120043012112345678

SN 12345678
PID 89
BIOTRONIK

SN 12345678
PID 89
BIOTRONIK



65bcedf

Acticor 7 HF-T QP DF4/IS4

SN 12345678

BIOTRONIK

Acticor 7 HF-T QP DF4/IS4

SN 12345678

BIOTRONIK

Acticor 7 HF-T QP DF4/IS4

SN 12345678

BIOTRONIK

Acticor 7 HF-T QP DF4/IS4

SN 12345678

BIOTRONIK

Acticor 7 HF-T QP DF4/IS4

REF 429522 PID 89

SN 12345678

BIOTRONIK

Acticor 7 HF-T QP DF4/IS4

REF 429522 PID 89

SN 12345678

BIOTRONIK

Acticor 7 HF-T QP DF4/IS4

REF 429522 SN 12345678 PID 89

BIOTRONIK

Acticor 7 HF-T QP DF4/IS4

REF 429522 SN 12345678 PID 89

BIOTRONIK

Acticor 7 HF-T QP DF4/IS4


REF 429522 SN 12345678 PID 89

BIOTRONIK

Acticor 7 HF-T QP


DF4/IS4 ProMRI Data Sheet

IS-1
IS4
DF4-LLHH



REF 429522 SN 12345678 2020-04-30 PID 89

Parameter	Values at Factory Release	Programmable Parameter Ranges
X-Ray Identification	+	-
Bradycardia Therapy Mode	OFF	AAI(R); V00; VVI(R); VVI-CLS; VDI(R); VDD(R); D00; DDI(R); DDD(R); DDD-ADI(R); DDD-CLS; OFF
Basic Rate; min-1	30...[5]...100...[10]...160	0.5...[0.25]...4...[0.5]...6; 7.5
Pulse Amplitude; V	0.4; 0.5...[0.25]...1.5	OFF; BURST; RAMP; HF BURST
Pulse Width; ms	OFF; BURST; RAMP	OFF; BURST; RAMP
A ATPs	OFF	OFF; DEFIBRILLATION
V ATPs	OFF	with/without Confirmation
Shock	OFF	



Manufacturer:
BIOTRONIK SE & Co. KG
Woermannkehre 1
12359 Berlin - Germany
Tel. +49 (0) 30 68905-0
Fax +49 (0) 30 6852804
sales@biotronik.com

Distributor:
BIOTRONIK, Inc.
6024 Jean Road
Lake Oswego, OR 97035-5369
(800) 547-0394 (24-hour)
(503) 635-9936 (Fax)

BIOTRONIK

441732

