

**BIOTRONIK**

**DF4/IS4**

**Intica 7 HF-T QP**

Use By 2017-10-31  
41°F / 113°F  
Stor.-Temp.

Serial No. SN 12345678

Cat.No. REF 404630



ID 101040354791421791717103112112345678



VDE-DDDRV  
Implantable Cardioverter Defibrillator

ProMRI®

**Intica 7 HF-T QP**

Prod.ID PID 89

Caution: Federal law (U.S.A.) restricts this device to sale by, or on the order of, a physician (or properly licenced practitioner)!



FCC ID: QR1TCHTNT2

426020

Manufacturer: BIOTRONIK SE & Co. KG  
Woermannkehrle 1  
12359 Berlin - Germany  
Tel. +49 (0) 30 68905-0  
Fax +49 (0) 30 6852804  
sales@biotronik.com

Distributor: BIOTRONIK, Inc.  
6024 Jean Road  
Lake Oswego, OR 97035-5369  
(800) 547-0394 (24-hour)  
(503) 635-9936 (Fax)

**BIOTRONIK**

**Intica 7 HF-T QP DF4/IS4**

REF 404630

Cat.No.

Use By 2017-10-31

SN 12345678

Serial No.

PID 89

Prod.ID



**BIOTRONIK**

Intica 7 HF-T QP DF4/IS4  
SN 12345678  
PID 89  
REF 404630  
Cat.No.


Intica 7 HF-T QP DF4/IS4  
SN 12345678  
PID 89  
REF 404630  
Cat.No.

Intica 7 HF-T QP DF4/IS4  
SN 12345678  
PID 89  
REF 404630  
Cat.No.

Intica 7 HF-T QP DF4/IS4  
SN 12345678  
PID 89  
REF 404630  
Cat.No.

**Intica 7 HF-T QP**

IS-1  
IS4-LLLL  
DF4-LLHH



DF4/IS4  
ProMRI®

Implantable Cardioverter Defibrillator

REF 404630

Cat.No.

Use By 2017-10-31


SN 12345678

Serial No.

PID 89

Prod.ID

STERILEEO  
Sterilized with ethylene oxide



For single use only, do not resterilize. Do not use if packaging is damaged. Refer to technical manual for instructions. ICD is shipped with all therapy disabled.

BIOTRONIK SE & Co. KG  
Woermannkehrle 1  
12359 Berlin - Germany

**BIOTRONIK**



ID 101040354791421791717103112112345678

SN 12345678  
PID 89  
REF 404630  
Cat.No.

SN 12345678  
PID 89  
REF 404630  
Cat.No.



c5034b2

Intica 7 HF-T QP DF4/IS4  
SN 12345678  
PID 89

**BIOTRONIK**

Intica 7 HF-T QP DF4/IS4  
SN 12345678  
PID 89

**BIOTRONIK**

Intica 7 HF-T QP DF4/IS4  
SN 12345678  
PID 89

**BIOTRONIK**

Intica 7 HF-T QP DF4/IS4  
SN 12345678  
PID 89

**BIOTRONIK**

**Intica 7 HF-T QP DF4/IS4**

REF 404630

PID 89

SN 12345678

**BIOTRONIK**

**Intica 7 HF-T QP DF4/IS4**

REF 404630

PID 89

SN 12345678

**BIOTRONIK**

**Intica 7 HF-T QP DF4/IS4**

REF 404630

SN 12345678

PID 89

**BIOTRONIK**

**Intica 7 HF-T QP DF4/IS4**

REF 404630

SN 12345678

PID 89

**BIOTRONIK**

**Intica 7 HF-T QP DF4/IS4**

REF 404630

SN 12345678

PID 89

**BIOTRONIK**

**Intica 7 HF-T QP**

**DF4/IS4 ProMRI® Data Sheet**

IS-1  
IS4-LLLL  
DF4-LLHH



REF 404630

SN 12345678

Use By 2017-10-31

PID 89

Implantable Cardioverter Defibrillator

Parameter	Values at Factory Release	Programmable Parameter Ranges
Bradycardia Therapy Mode	OFF	AAI(R); V00; VVI(R); VVI-CLS; VDDI(R); VDD(R); D00; DDI(R); DDD(R); DDD-ADI(R); DDD-CLS; OFF
Basic Rate; min-1		30...[5]...100...[10]...160
Pulse Amplitude; V		0.5...[0.25]...4...[0.5]...6; 7.5
Pulse Width; ms		0.4; 0.5...[0.25]...1.5
A ATPs	OFF	OFF; BURST; RAMP; HF BURST
V ATPs	OFF	OFF; BURST; RAMP
Shock	OFF	OFF; DEFIBRILLATION with/without Confirmation
X-Ray Identification	NK	-

For single use only. Refer to Technical Manual for Instructions.

Distributor: BIOTRONIK, Inc.  
6024 Jean Road  
Lake Oswego, OR 97035-5369  
(800) 547-0394 (24-hour)  
(503) 635-9936 (Fax)

BIOTRONIK SE & Co. KG  
Woermannkehrle 1  
12359 Berlin - Germany  
Tel. +49 (0) 30 68905-0  
Fax +49 (0) 30 6852804  
sales@biotronik.com



426020

**BIOTRONIK**