

# APPENDIX A – SAMPLE LABEL

## Labelling Requirements

The sample label shown shall be *permanently affixed* at a conspicuous location on the device and be readily visible to the user at the time of purchase.

Distributed by:  
**stryker®**  
**ENDOSCOPY**  
Model : AGM19TK  
REF : 240-030-900  
19" Strykerviison 1  
Flat Panel Monitor  
  
DC 12V --- 6.67A  
S/N :  
QUANTITY : 1  
  
5900 OPTICAL COURT  
SAN JOSE, CA 95138 U.S.A.

⚠

AC/DC Adapter

Manufacturer : Ault Inc.

or Stryker Model : 240-030-901

⚠

CAUTION - See accompanying documents.

⚠

CAUTION - ELECTRIC SHOCK

"To avoid electrical shock, do not open the cabinet.

Refer servicing to qualified personnel only."

CLASSIFIED  
C  
UL  
US  
51LJ  
Medical Equipment  
E215822

UL 2601-1

CAN/CSA C22.2

NO 501.1

FCC ID : QVXAGM19TK

This device is compliant with Part 15 of the FCC Rules. Operation is subject to the following two conditions: (1) this device may not cause harmful interference, and (2) this device must accept any interference received, including interference that may cause.

Advan International Corp.

● **FCC ID Location of EUT**

ADVAN International Corp.  
FCC ID:QVXAGM19TK

Page 21 of 40