

EDAN Blood Gas and Chemistry Analyzer

Model: i15



YYYY-MM

SN



CCC

Anti-electro-shock Type: Class I

100V-240V~ 50Hz/60Hz 1.2A-0.5A

FCC ID:SMQI15EDAN

EC REP

Shanghai international Holding Corp.GmbH(Europe)
Eiffestrasse 80,D-20537 Hamburg Germany

Tel:+49-40-2513175 Fax:+49-40-255726



IVD



Rx only (U.S.)



Edan Instruments, Inc.

3/F - B,Nanshan Medical Equipments Park,Nanghai Rd 1019#,

Shekou,Nanshan Shenzhen,518067 P.R. China

Tel:+86-755-26898326 Fax:+86-755-26898330

[http:// www.edan.com.cn](http://www.edan.com.cn)

Label:

FCC ID : SMQI15EDAN

This device complies with Part 15 of the FCC Rules. Operation is subject to the following two conditions: (1) this device may not cause harmful interference, and (2) this device must accept any interference received, including interference that may cause undesired operation.

Label Location:

