


41.4 mm



This device complies with Part 15 of the FCC Rules. Operation is subject to the following two conditions:
(1) This device may not cause harmful interference, and
(2) this device must accept any interference received, including interference that may cause undesired operation.

EDAN Electrocardiograph

Model: iSE-301

YYYY-MM

CE 0123

SN



CCCC

Rx only

19V 2A

14.8V Rechargeable Li-ion Battery

FCC ID: SMQSE301EDAN

EC REP Shanghai International Holding Corp.GmbH(Europe)
Eiffestrasse 80,D-20537 Hamburg Germany
Tel:+49-40-2513175 Fax:+49-40-255726



Edan Instruments, Inc.
3/F - B,Nanshan Medical Equipments Park,Nanghai
Rd 1019#, Shekou,Nanshan Shenzhen,518067 P.R. China
Tel:+86-755-26898326 Fax:+86-755-26898330
[http:// www.edan.com.cn](http://www.edan.com.cn)

103.6 mm