

Angel Medical Systems, Inc.

Guardian[®] IMD/EXD System for monitoring electrograms (EGMs), 402 MHz – 405 MHz

FCC ID: THL- IMDAG101

Photo Four – External

Implantable Medical Device (IMD) Printed Art Work of IMD Implant Sterile Packaging with FCC ID Number ... THL-IMDAG101

Note: The IMD is a Sterile Product, implanted and not visible during operation.

AngelMed Guardian[®] IMD Model # AG101



Implantable Medical Device (IMD) for detecting and warning a patient of a variety of changes in his or her electrogram.

For use ONLY with Angel Medical Systems, Inc. Model LA-001 Lead Adapter and St. Jude Model 1488T Bipolar Cardiac Lead

Contents: One IMD, One Torque Wrench

The contents (STERILE) of this inner package have been sterilized by ethylene-oxide gas. Sterility cannot be guaranteed if sterile package is damaged or opened.



STERILE EO



YYYY-MM



YYYY-MM



SN

XXXXXX



-10°C/14° F and +55°C/131°F

Caution: Investigational Device. Limited by federal (or United States) law to investigational use.

Manufactured for: Angel Medical Systems, Inc.
1 Sheila Drive
Tinton Falls, NJ 07724 USA
Phone: 732-212-1888 (USA)



THL-IMDAG101

(Hi-Tronics Document No.)