



## **Test Setup Photos-FCC**

|                                    |                                                                                                                                                          |
|------------------------------------|----------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>Applicant:</b>                  | Beijing Choice Electronic Technology Co., Ltd.                                                                                                           |
| <b>Address of Applicant:</b>       | 2nd Floor, 3rd Floor and Room 410-412 4th Floor, No. 2 Building, No. 9 Shuangyuan Road, Shijingshan District, 100041 Beijing, PEOPLE'S REPUBLIC OF CHINA |
| <b>Equipment Under Test (EUT):</b> |                                                                                                                                                          |
| <b>Product Name:</b>               | Wrist Pulse Oximeter                                                                                                                                     |
| <b>Model No.:</b>                  | MD300W628                                                                                                                                                |




## **1 Contents**

|                                              | Page     |
|----------------------------------------------|----------|
| <b>1 CONTENTS.....</b>                       | <b>2</b> |
| <b>2 EQUIPMENT UNDER TEST PICTURES .....</b> | <b>3</b> |
| 2.1 Radiated Emission Test Setup .....       | 3        |
| 2.2 RF Conducted Test Setup .....            | 4        |

## 2 EQUIPMENT UNDER TEST PICTURES

### 2.1 Radiated Emission Test Setup

30MHz to 1GHz




Above 1GHz



Zoom:



## 2.2 RF Conducted Test Setup



**End**